




Optimize your distillation

With the BUCHI Multivapor™



Remember: for evaporation at 30 °C (vapor temperature) follow the guideline:

1. instrument setting: 55 °C
2. cooling medium: ≤ 10 °C, ≥ 500 mL/min
3. vacuum setting according to the Solvent List column 30 °C

Solvent	Formula	Vacuum [mbar] for bp at:	
		30 °C	50 °C
Acetic acid	C ₂ H ₄ O ₂	26	72
Acetone	C ₃ H ₆ O	370	815
Acetonitrile	C ₂ H ₃ N	134	315
<i>n</i> -Amyl alcohol, <i>n</i> -Pentanol	C ₅ H ₁₂ O	6	20
<i>n</i> -Butanol	C ₄ H ₁₀ O	14	44
<i>tert</i> -Butanol	C ₄ H ₁₀ O	78	231
Chlorobenzene	C ₆ H ₅ Cl	22	56
Chloroform	CHCl ₃	306	665
Cyclohexane	C ₆ H ₁₂	154	347
1,2-Dichloroethane	C ₂ H ₄ Cl ₂	137	315
Dichloromethane	CH ₂ Cl ₂	685	atm. press.
Diethylether	C ₄ H ₁₀ O	838	atm. press.
<i>trans</i> -1,2-Dichloroethylene	C ₂ H ₂ Cl ₂	317	705
Diisopropylether	C ₆ H ₁₄ O	251	545
Dioxane	C ₄ H ₈ O ₂	68	165
Dimethylformamide (DMF)	C ₃ H ₇ NO	6	17
Ethanol	C ₂ H ₆ O	97	276
Ethylacetate	C ₄ H ₈ O ₂	153	366
Heptane	C ₇ H ₁₆	77	183
Hexane	C ₆ H ₁₄	241	525
Isopropylalcohol	C ₃ H ₈ O	78	231
Isoamylalcohol	C ₅ H ₁₂ O	9	29
Methyl <i>tert</i> -butyl ether (MTBE)	C ₅ H ₁₂ O	413	835
Methyl ethyl ketone (MEK)	C ₄ H ₈ O	160	359
Methanol	CH ₄ O	236	607
Pentane	C ₅ H ₁₂	819	atm. press.
<i>n</i> -Propanol	C ₃ H ₈ O	37	115
Pentachloroethane	C ₂ HCl ₅	8	21
Tetrachloromethane	CCl ₄	179	398
Tetrahydrofuran (THF)	C ₄ H ₈ O	234	539
Toluene	C ₇ H ₈	48	118
Trichloroethylene	C ₂ HCl ₃	119	275
Water	H ₂ O	42	120
Xylene	C ₈ H ₁₀	15	40

